

# Orders API

The Orders endpoint represents a resource containing approved requests. Requests specify a set of samples approved for processing in Lab Operations. The endpoint parameters are applied to the operations order retrieval, creation, cancellation, reset, completion and update.

Relevant endpoints are described in each topic heading, required and optional parameters, typical requests and responses (in JSON format) are described following. Optional parameters and objects are listed in a reference.

- Order
  - Order - standard parameters
    - order Example - JSON format
  - Query Orders (objects)
    - getOrder Response Example - JSON Format
  - Create Order
    - Create Order Example - JSON Format
  - Create Orders - batch
  - Cancel Orders
  - Reset Orders
  - Complete Orders
  - Update Order
  - Order Objects
    - attachment
    - comments
    - recipients
    - format
    - lineItems
    - intermediateFormat
    - comment
    - miracleRun
  - Nested Attributes
    - wellVectors
    - sourceInventoryItems
    - orderFragments
    - selectedPath

## Order

- Methods:
  - GET: retrieve a list of orders
  - POST: create (many)
  - PUT: update, cancel, complete, reset
  - DELETE: remove
- Headers:
  - application/json
  - Pagination
    - skip (number): parameter that defines the number of entries that are ignored before returning a response.
    - take (number): parameter that defines the number of entries returned in a single response.
    - count (boolean): true or false boolean.

## Order - standard parameters

A *well formed* order in Lab Ops specifies the parameters described in the following table.

<sup>1</sup> Orders id value is generated automatically by Lab Ops for orders in the database.

#	Parameter	Parent	Type	Example
1	id <sup>1</sup>	-	UUID	01z2GaZzKuS3s5cX6npsKD
2	name	-	string	test 1
3	description	-	string	...
4	tags	-	array	[ ]
5	recipients	-	<b>object</b>	success@simulator. amazonsees.com

6	recipients address	-	string	Spinnereistraße 7, 95445 Bayreuth
7	format	-	<b>object</b>	vial liquid
8	name	format	string	384 - IC50 - 8 Points - 3.162 DF
9	targetContainerType	format	string	GR384
10	materialState	format	string	Liquid
11	dilutionSeries	format	numbers (array)	[3.162, 3.162, 3.162 ...]
12	wellVectors	<b>format</b>	<b>object</b>	?
13	type	<i>wellVectors</i>	string	Sample
14	instanceNumber	<i>wellVectors</i>	number	0
15	start	<i>wellVectors</i>	<b>object</b>	-
16	rows	<i>start</i>	number	0
17	columns	<i>start</i>	number	0
18	wells	<i>wellVectors</i>	number	1
19	direction	<i>wellVectors</i>	string	East
20	spacing	<i>wellVectors</i>	number	0
21	index	<i>wellVectors</i>	<b>number</b>	0
22	rows	format	number	16
23	columns	format	number	24
24	lineItems	-	<b>object</b>	-
25	amount	<b>lineItems</b>	number	60
26	topConcentration	<b>lineItems</b>	number	10
27	status	<b>lineItems</b>	string	pending
28	sourceType	<b>lineItems</b>	string	sample
29	sourceInventoryItem	<b>lineItems</b>	<b>object</b>	-
30	containerBarcode	sourceInventoryItem	string	Test tube
31	position	<b>sourceInventoryItem</b>	<b>object</b>	-
32	row	position	array	0
33	column	position	array	0
34	externalReferences	lineItems	array	[ ]
35	orderFragments	-	array	[ ]
36	status	-	date	created
37	isLineSequenceRelevant	-	string	false
38	controlVolume	-	number	1
39	numberOfCopies	-	number	1

### order Example - JSON format

```
{
  "name": "Late binding test 1",
  "description": "Order to test late binding",
  "tags": [],
  "recipients": ["success@simulator.amazonses.com"],
  "recipientAddress": "Cypress",
```

```

"format": {
  "name": "Solubilization Test Format",
  "targetContainerType": "1.4ml Tube",
  "materialState": "Liquid",
  "dilutionSeries": [],
  "wellVectors": [
    {
      "type": "Sample",
      "instanceNumber": 1,
      "start": {
        "row": 0,
        "column": 0
      },
      "wells": 0,
      "direction": "East",
      "spacing": 0,
      "index": 0
    }
  ],
  "rows": 1,
  "columns": 1
},
"lineItems": [
  {
    "amount": 60,
    "topConcentration": 10,
    "status": "Pending",
    "sourceType": "Sample",
    "sourceInventoryItem": {
      "containerBarcode": "LateBindingTestTube",
      "position": {
        "row": 0,
        "column": 0
      }
    },
    "externalReferences": []
  }
],
"orderFragments": [],
"status": "Created",
"isLineItemSequenceRelevant": false,
"controlVolume": 1,
"numberOfCopies": 1
}

```

#### Query Orders (objects)

- Method:
  - GET
- lab-operations/{lab}/orders2

Query parameters for the `orders2` endpoint specifies orders retrieval in the Lab Operations database and are described in the following table.

<sup>2</sup> Denotes attributes that apply to every query object except for `id` and `type`.

#	Parameter	Parent	Type	Example
1	<code>id</code>	-	<b>object</b>	4
2	<code>mode</code>	<code>id</code>	<b>boolean</b>	default
3	<code>in</code>	<code>id</code>	string (array)	-
4	<code>notIn</code>	<code>id</code>	string (array)	-
5	<code>name</code>	-	<b>object</b>	-
6	<code>equals<sup>2</sup></code>	<code>name</code>	string	-
7	<code>mode<sup>2</sup></code>	<code>name</code>	<b>boolean</b>	default
8	<code>in<sup>2</sup></code>	<code>name</code>	string (array)	-
9	<code>notIn<sup>2</sup></code>	<code>name</code>	string (array)	-
10	<code>contains<sup>2</sup></code>	<code>name</code>	string	-
11	<code>startsWith<sup>2</sup></code>	<code>name</code>	string	-
12	<code>endsWith<sup>2</sup></code>	<code>name</code>	string	7
13	<code>format</code>	-	<b>object</b>	Liquid
14	<code>status</code>	-	<b>object</b>	inactive
15	<code>type</code>	-	<b>object</b>	-
16	<code>in</code>	<code>type</code>	string (array)	-
17	<code>value</code>	-	<b>string</b>	-
18	<code>fields</code>	-	string (array)	name, status
19	<code>message</code>	-	string	-
20	<code>tags</code>	-	<b>object</b>	-
21	<code>lab</code>	-	string (path)	bt

### getOrder Response Example - JSON Format

The following code snippet demonstrates the response returned when a query is submitted with some values left empty.

▼ [getOrder Response JSON Format](#)

```
{
  "orders": [
    {
      "tags": [],
      "priority": "Low",
      "numberOfCopies": 1,
      "id": "5788915c-af39-46a1-90ea-5c6d4a1cb61f",
      "name": "Star-Echo + ACI",
      "description": "Star Echo Order with ACI",
      "createdBy": "System",
      "recipients": [
        "success@simulator.amazonses.com",
        "john@smith.com",
        "james@hetfield.org"
      ],
    },
  ],
}
```

```
"recipientAddress": "Spinnereistraße 7, 95445 Bayreuth",
"receivedAt": "2022-07-07T10:23:16.296Z",
"status": "Ready",
"hasStartedJobs": false,
"allLineItemsHaveDeliverable": false
},
{
  "tags": [],
  "priority": "High",
  "numberOfCopies": 1,
  "id": "c5b588b1-37a9-46b8-a09d-1763f08849ff",
  "name": "Serialize from preintermediate",
  "description": "this one models a serialization from a
preintermediate (Star-Star)",
  "createdBy": "System",
  "recipients": [
    "success@simulator.amazonses.com"
  ],
  "recipientAddress": "Spinnereistraße 7, 95445 Bayreuth",
  "receivedAt": "2022-07-07T10:23:14.118Z",
  "status": "Ready",
  "hasStartedJobs": false,
  "allLineItemsHaveDeliverable": false
},
{
  "tags": [],
  "priority": "High",
  "numberOfCopies": 1,
  "id": "09a00e19-ffca-4bbb-9483-f043f2db7215",
  "name": "Test Order 3",
  "description": "Transfer SourceTube to Plate",
  "createdBy": "System",
  "recipients": [
    "success@simulator.amazonses.com"
  ],
  "recipientAddress": "Spinnereistraße 7, 95445 Bayreuth",
  "receivedAt": "2022-07-07T10:23:14.806Z",
  "status": "Ready",
  "hasStartedJobs": false,
  "allLineItemsHaveDeliverable": false
},
{
  "tags": [],
  "priority": "Medium",
  "numberOfCopies": 1,
  "id": "aa114394-190b-4e29-ae4d-2b9c2a3812da",
  "name": "New Order 2",
  "description": "Sources on Plates",
  "createdBy": "System",
  "recipients": [
```

```
    "success@simulator.amazonses.com"
  ],
  "recipientAddress": "Spinnereistraße 7, 95445 Bayreuth",
  "receivedAt": "2022-07-07T10:23:15.167Z",
  "status": "Ready",
  "hasStartedJobs": false,
  "allLineItemsHaveDeliverable": false
},
{
  "tags": [],
  "priority": "Low",
  "numberOfCopies": 1,
  "id": "16d8c109-573a-4b30-9701-36ddd2dbf18a",
  "name": "Solubilization + Star Only",
  "description": "Solubilization and star only order",
  "createdBy": "System",
  "recipients": [
    "success@simulator.amazonses.com"
  ],
  "recipientAddress": "Spinnereistraße 7, 95445 Bayreuth",
  "receivedAt": "2022-07-07T10:23:17.291Z",
  "status": "Ready",
  "hasStartedJobs": false,
  "allLineItemsHaveDeliverable": false
},
{
  "tags": [],
  "priority": "Low",
  "numberOfCopies": 1,
  "id": "cbfe5689-b45c-4f3a-b155-99dc1126bb5b",
  "name": "ValidateSourceLiquidOrder",
  "description": "validate source order",
  "createdBy": "Jason",
  "recipients": [
    "success@simulator.amazonses.com"
  ],
  "recipientAddress": "Spinnereistraße 7, 95445 Bayreuth",
  "receivedAt": "2022-09-20T08:24:02.711Z",
  "status": "Released",
  "hasStartedJobs": false,
  "allLineItemsHaveDeliverable": false
},
{
  "tags": [],
  "priority": "High",
  "numberOfCopies": 1,
  "id": "559f3024-3da1-4828-b489-50de77492318",
  "name": "Neat order",
  "description": "neat order",
  "createdBy": "System",
```

```
"recipients": [
  "success@simulator.amazonses.com"
],
"recipientAddress": "Spinnereistraße 7, 95445 Bayreuth",
"receivedAt": "2022-07-07T10:23:17.850Z",
"status": "Ready",
"hasStartedJobs": false,
"allLineItemsHaveDeliverable": false
},
{
  "tags": [],
  "priority": "Medium",
  "numberOfCopies": 1,
  "id": "belc20e4-bcda-47f2-85bb-3876caca02ea",
  "name": "Big Order 1",
  "description": "Transfer from rack to plate",
  "createdBy": "System",
  "recipients": [
    "success@simulator.amazonses.com"
  ],
  "recipientAddress": "Spinnereistraße 7, 95445 Bayreuth",
  "receivedAt": "2022-07-07T10:23:16.602Z",
  "status": "Released",
  "hasStartedJobs": false,
  "allLineItemsHaveDeliverable": false
},
{
  "tags": [],
  "priority": "Low",
  "numberOfCopies": 1,
  "id": "1de156d2-23ee-46bc-ae7f-f55ee51cd57c",
  "name": "Star Only + Reference + Dilution Series",
  "description": "Order with references and dilution series",
  "createdBy": "System",
  "recipients": [
    "success@simulator.amazonses.com"
  ],
  "recipientAddress": "Spinnereistraße 7, 95445 Bayreuth",
  "receivedAt": "2022-07-07T10:23:15.574Z",
  "status": "Ready",
  "hasStartedJobs": false,
  "allLineItemsHaveDeliverable": false
},
{
  "tags": [],
  "priority": "High",
  "numberOfCopies": 1,
  "id": "7e2faf5d-9c9b-4c64-b53f-2e5a68d03abf",
  "name": "Test Order 4",
  "description": "Transfer SourceTubes from ManStore to Plate",
```

```
"createdBy": "System",
"recipients": [
  "success@simulator.amazonses.com"
],
"recipientAddress": "Spinnereistraße 7, 95445 Bayreuth",
"receivedAt": "2022-07-07T10:23:14.958Z",
"status": "Ready",
"hasStartedJobs": false,
"allLineItemsHaveDeliverable": false
},
{
  "tags": [],
  "priority": "Low",
  "numberOfCopies": 1,
  "id": "763af6b5-e826-4517-9722-722faf95ale3",
  "name": "Star-Bravo + Reference",
  "description": "Order with references",
  "createdBy": "System",
  "recipients": [
    "success@simulator.amazonses.com"
  ],
  "recipientAddress": "Spinnereistraße 7, 95445 Bayreuth",
  "receivedAt": "2022-07-07T10:23:15.783Z",
  "status": "Ready",
  "hasStartedJobs": false,
  "allLineItemsHaveDeliverable": false
},
{
  "tags": [],
  "priority": "Low",
  "numberOfCopies": 1,
  "id": "f9c84719-1390-4ee8-b7c9-9471ca102d5b",
  "name": "New Order 3",
  "description": "Sources on plates 2",
  "createdBy": "System",
  "recipients": [
    "success@simulator.amazonses.com"
  ],
  "recipientAddress": "Spinnereistraße 7, 95445 Bayreuth",
  "receivedAt": "2022-07-07T10:23:15.392Z",
  "status": "Ready",
  "hasStartedJobs": false,
  "allLineItemsHaveDeliverable": false
},
{
  "tags": [],
  "priority": "High",
  "numberOfCopies": 1,
  "id": "08626ce9-6f5c-40a8-9b50-97986198456a",
  "name": "Test Order",
```



```
"description": "Transfer SourceTube to Plate",
"createdBy": "System",
"recipients": [
  "success@simulator.amazonses.com"
],
"recipientAddress": "Spinnereistraße 7, 95445 Bayreuth",
"receivedAt": "2022-07-07T10:23:14.421Z",
"status": "No path",
"hasStartedJobs": false,
"allLineItemsHaveDeliverable": false
},
{
  "tags": [],
  "priority": "High",
  "numberOfCopies": 1,
  "id": "cb25deb0-8f0a-4be0-aedb-492741b02bd6",
  "name": "Test Order 2",
  "description": "Transfer SourceTube to Plate",
  "createdBy": "System",
  "recipients": [
    "success@simulator.amazonses.com"
  ],
  "recipientAddress": "Spinnereistraße 7, 95445 Bayreuth",
  "receivedAt": "2022-07-07T10:23:14.617Z",
  "status": "Ready",
  "hasStartedJobs": false,
  "allLineItemsHaveDeliverable": false
},
{
  "tags": [],
  "priority": "Low",
  "numberOfCopies": 1,
  "id": "def1c6e5-5b4e-459d-a4ab-e8a27800e2d9",
  "name": "Star Only + different volumes and concentrations",
  "description": "Order with references and dilution series",
  "createdBy": "System",
  "recipients": [
    "success@simulator.amazonses.com"
  ],
  "recipientAddress": "Spinnereistraße 7, 95445 Bayreuth",
  "receivedAt": "2022-07-07T10:23:16.126Z",
  "status": "Ready",
  "hasStartedJobs": false,
  "allLineItemsHaveDeliverable": false
},
{
  "tags": [],
  "priority": "High",
  "numberOfCopies": 1,
  "id": "2764a14a-fbaa-4108-a8a7-b32ef73661ac",
```

```
    "name": "Echo order (no ACI)",
    "description": "Dilute from intermediate order that doesn't
require an ACI",
    "createdBy": "System",
    "recipients": [
      "success@simulator.amazonses.com"
    ],
    "recipientAddress": "Spinnereistraße 7, 95445 Bayreuth",
    "receivedAt": "2022-07-07T10:23:17.557Z",
    "status": "Ready",
    "hasStartedJobs": false,
    "allLineItemsHaveDeliverable": false
  },
  {
    "tags": [],
    "priority": "Low",
    "numberOfCopies": 1,
    "id": "cb7c43f9-aa97-44b3-80f4-7d4212568e90",
    "name": "Split Neat 0014",
    "description": "order for split neat job",
    "createdBy": "System",
    "recipients": [
      "success@simulator.amazonses.com",
      "jimmy@hendrix.org",
      "alpha@centauri.net"
    ],
    "recipientAddress": "Spinnereistraße 7, 95445 Bayreuth",
    "receivedAt": "2022-07-07T10:23:15.209Z",
    "status": "Released",
    "hasStartedJobs": false,
    "allLineItemsHaveDeliverable": false
  },
  {
    "tags": [],
    "priority": "High",
    "numberOfCopies": 1,
    "id": "de7c0537-ba3b-463d-b70e-b57495cd7280",
    "name": "No Path Order",
    "description": "To be tested with miracle run",
    "createdBy": "System",
    "recipients": [
      "success@simulator.amazonses.com"
    ],
    "recipientAddress": "Spinnereistraße 7, 95445 Bayreuth",
    "receivedAt": "2022-07-07T10:23:17.421Z",
    "status": "No path",
    "hasStartedJobs": false,
    "allLineItemsHaveDeliverable": false
  },
  {
```

```

    "tags": [],
    "priority": "Low",
    "numberOfCopies": 1,
    "id": "4c824bf2-27b8-40a8-82da-6c92debd4071",
    "name": "Star-Echo + Reference",
    "description": "Star Echo Order with references",
    "createdBy": "System",
    "recipients": [
      "success@simulator.amazonses.com"
    ],
    "recipientAddress": "Spinnereistraße 7, 95445 Bayreuth",
    "receivedAt": "2022-07-07T10:23:15.940Z",
    "status": "Ready",
    "hasStartedJobs": false,
    "allLineItemsHaveDeliverable": false
  },
  {
    "tags": [],
    "priority": "High",
    "numberOfCopies": 1,
    "id": "aa20d3bb-9f2b-4fc7-bc58-af20eb64ea4a",
    "name": "Serialization order (wrapped dilution series)",
    "description": "order that has a dilution series wrapped around
two rows",
    "createdBy": "System",
    "recipients": [
      "success@simulator.amazonses.com"
    ],
    "recipientAddress": "Spinnereistraße 7, 95445 Bayreuth",
    "receivedAt": "2022-07-07T10:23:17.679Z",
    "status": "Ready",
    "hasStartedJobs": false,
    "allLineItemsHaveDeliverable": false
  }
]
}

```

## Create Order

- Method
  - POST
- `lab-operations/{lab}/orders`

The `orders` endpoint allows you to create a structured request that defines a material sample to be processed by Lab Operations.

Objects can be specified individually or batch. Object parameters that can be applied to an order specification are described in the following table.

#	Parameter	Parent	Type/	Example
1	<code>id</code>	-	string	<code>01z2GaZzKuS3s5cX6npsKD</code>
2	<code>name</code>	-	string	<code>test 1</code>

3	description	-	string	an order
4	createdBy	-	string	userID?
5	createdAt	-	date	2022-11-11T16:14:43.692Z
6	updatedBy	-	string	-
7	updatedAt	-	date	2022-11-11T16:14:43.692Z
8	attachments	-	object	-
9	tags	-	array	november
10	comments	-	object (array)	[ ]
11	recipients	-	object (array)	success@simulator.amazonses.com
12	recipients address	-	string	Spinnereistraße 7, 95445 Bayreuth
13	receivedAt	-	string	2022-11-11T16:14:43.692Z
14	format	-	object	-
15	lineItems	-	object	-
16	intermediateFormat	lineItems	object	-
17	sourceInventoryItem	lineItems	object	-
18	orderFragments	-	object (array)	-
19	status	-	date	created
20	type	-	object	Fulfillment
21	isNotExecutable	-	boolean	true
22	priority	-	object	-
23	numberOfCopies	-	number	1
24	isLineSequenceRelevant	-	string	false
25	controlVolume	-	number	1
26	miracleRun	-	object	-
26	message	-	string	requested volume has...
27	hasStartedJobs	-	boolean	true
28	allLineItemsHaveDeliverable	-	boolean	true
29	numberOfLineItems	-	number	0
30	numberOfDeliverables	-	number	1
31	dueDate	-	date	2022-11-17T12:25:55.388Z

## Create Order Example - JSON Format

▼ [createOrders JSON Format](#)

```
{
  "id": "string",
  "name": "string",
  "description": "string",
  "createdBy": "string",
  "createdAt": "2022-11-17T12:25:55.387Z",
  "updatedBy": "string",
  "updatedAt": "2022-11-17T12:25:55.387Z",
  "attachments": [
    {
```

```
    "id": "string",
    "name": "string",
    "previewEnabled": true,
    "role": {},
    "size": 0,
    "mimeType": "string",
    "payload": {}
  }
],
"tags": [
  "string"
],
"comments": [
  {
    "id": "string",
    "value": "string",
    "timestamp": "2022-11-17T12:25:55.388Z",
    "userDisplayName": "string"
  }
],
"recipients": [
  "string"
],
"recipientAddress": "string",
"receivedAt": "2022-11-17T12:25:55.388Z",
"format": {
  "id": "string",
  "name": "string",
  "targetContainerType": "string",
  "targetBarcode": "string",
  "targetAmount": 0,
  "materialState": {},
  "wellVectors": [
    {
      "id": "string",
      "instanceNumber": 0,
      "type": {},
      "index": 0,
      "start": {
        "row": 0,
        "column": 0
      },
      "wells": 0,
      "direction": {},
      "spacing": 0
    }
  ],
  "dilutionSeries": [
    null
  ],

```

```
"rows": 0,
"columns": 0
},
"lineItems": [
  {
    "id": "string",
    "status": "string",
    "amount": 0,
    "topConcentration": 0,
    "recipients": [
      "string"
    ],
    "recipientAddress": "string",
    "intermediateFormat": {
      "id": "string",
      "name": "string",
      "targetContainerType": "string",
      "targetBarcode": "string",
      "targetAmount": 0,
      "materialState": {},
      "wellVectors": [
        {
          "id": "string",
          "instanceNumber": 0,
          "type": {},
          "index": 0,
          "start": {
            "row": 0,
            "column": 0
          },
          "wells": 0,
          "direction": {},
          "spacing": 0
        }
      ],
      "dilutionSeries": [
        null
      ],
      "rows": 0,
      "columns": 0
    },
    "sourceInventoryItem": {
      "sampleId": "string",
      "containerBarcode": "string",
      "location": {
        "id": "string",
        "displayName": "string",
        "type": {},
        "parentId": "string",
        "barcode": "string",
```

```
    "capacity": 0,
    "layout": [
      "string"
    ],
    "orientation": {},
    "countingDirection": {},
    "loadingSequence": {},
    "children": [
      "string"
    ],
    "locationTree": [
      "string"
    ],
    "hidden": true,
    "allowHazardous": true,
    "allowControlled": true
  },
  "position": {
    "row": 0,
    "column": 0
  }
},
"comment": "string",
"layoutIndex": 0,
"deliveryIndex": 0,
"externalReferences": [
  "string"
],
"sourceType": {},
"deliverableSampleId": "string"
}
],
"orderFragments": [
  {
    "id": "string",
    "order": {
      "id": "string",
      "name": "string",
      "type": {},
      "description": "string",
      "recipients": [
        "string"
      ],
      "priority": {},
      "receivedAt": "2022-11-17T12:25:55.388Z",
      "recipientAddress": "string",
      "numberOfCopies": {},
      "isLineItemSequenceRelevant": true
    },
  },
  "lineItems": [
```

```
{
  "id": "string",
  "status": "string",
  "amount": 0,
  "topConcentration": 0,
  "recipients": [
    "string"
  ],
  "recipientAddress": "string",
  "intermediateFormat": {
    "id": "string",
    "name": "string",
    "targetContainerType": "string",
    "targetBarcode": "string",
    "targetAmount": 0,
    "materialState": {},
    "wellVectors": [
      {
        "id": "string",
        "instanceNumber": 0,
        "type": {},
        "index": 0,
        "start": {
          "row": 0,
          "column": 0
        },
        "wells": 0,
        "direction": {},
        "spacing": 0
      }
    ],
    "dilutionSeries": [
      null
    ],
    "rows": 0,
    "columns": 0
  },
  "sourceInventoryItem": {
    "sampleId": "string",
    "containerBarcode": "string",
    "location": {
      "id": "string",
      "displayName": "string",
      "type": {},
      "parentId": "string",
      "barcode": "string",
      "capacity": 0,
      "layout": [
        "string"
      ]
    }
  ],
}
```



```
        "orientation": {},
        "countingDirection": {},
        "loadingSequence": {},
        "children": [
            "string"
        ],
        "locationTree": [
            "string"
        ],
        "hidden": true,
        "allowHazardous": true,
        "allowControlled": true
    },
    "position": {
        "row": 0,
        "column": 0
    }
},
"comment": "string",
"layoutIndex": 0,
"deliveryIndex": 0,
"externalReferences": [
    "string"
],
"sourceType": {},
"deliverableSampleId": "string"
}
],
"selectedPath": {
    "consumed": 0,
    "effort": 0,
    "runs": [
        "string"
    ],
    "strategy": {}
}
},
"status": "string",
"type": "Fulfillment",
"isNotExecutable": true,
"priority": {},
"numberOfCopies": {},
"isLineItemSequenceRelevant": true,
"controlVolume": 0,
"miracleRun": {
    "id": "string",
    "orderId": "string",
    "format": {
        "id": "string",
```

```

"name": "string",
"targetContainerType": "string",
"targetBarcode": "string",
"targetAmount": 0,
"materialState": {},
"wellVectors": [
  {
    "id": "string",
    "instanceNumber": 0,
    "type": {},
    "index": 0,
    "start": {
      "row": 0,
      "column": 0
    },
    "wells": 0,
    "direction": {},
    "spacing": 0
  }
],
"dilutionSeries": [
  null
],
"rows": 0,
"columns": 0
},
"amount": 0,
"topConcentration": 0,
"equipments": [
  "string"
],
"instructions": "string"
},
"message": "string",
"hasStartedJobs": true,
"allLineItemsHaveDeliverable": true,
"numberOfLineItems": 0,
"numberOfDeliverables": 0,
"dueDate": "2022-11-17T12:25:55.388Z"
}

```

### Create Orders - batch

- Method
  - POST
- /lab-operations/{lab}/orders/create-many

The `create-many` endpoint allows you to *batch* create structured requests that define material samples to be processed by Lab Operations.

### Cancel Orders

- Method
  - DELETE (PUT?)
- /lab-operations/{lab}/orders/cancelation

The cancelation endpoint allows you to cancel (delete) order by id in the Lab Operations database.

### Reset Orders

- Method
  - PUT
- /lab-operations/{lab}/orders/reset

The reset endpoint allows you to reset an order by id to its original status in the Lab Operations database.

### Complete Orders

- Method
  - PUT

/lab-operations/{lab}/orders/completion

The complete endpoint allows you to batch complete orders by id in bulk in the Lab Operations database.

### Update Order

- Method
  - PUT

/lab-operations/{lab}/orders/{id}

The {id} endpoint allows you to update the order directly in the Lab Operations database.

## Order Objects

### attachment

The attachment object is a container for images and document components of the order. These items can be included with the order from the request system or directly uploaded from a network or local machine.

#	Parameter	Parent	Type/ Reference	Example
1	attachments	-	object (array)	-
2	id	attachments	string	-
3	name	attachments	string	-
4	previewEnabled	attachments	boolean	true
5	role	attachments	object	-
6	size	attachments	number	0
7	mimeType	attachments	string	-
8	payload	attachments	object	-

### comments

The comments object is a container for a string displayed on the interface intended as instructions or information for the operator regarding the order. This item can be included with the order from the request system or directly uploaded from a network or local machine.

#	Parameter	Parent	Type	Example
1	comments	-	object (array)	[ ]
2	id	comments	string	-

3	value	comments	string	-
4	timestamp	comments	date	2022-11-11T16:14:43.692Z
5	userDisplayName	comments	string	-

## recipients

The `recipients` parameter is a string that contains the email address or name of a person or entity for whom the complete processed order is transferred. The `recipientAddress` is typically the physical address applied to the shipping label. The recipient is also informed in some configurations when certain process steps are complete.

#	Parameter	Parent	Type	Example
1	<code>recipients</code>	-	string (array)	john.doe@assayLab.com
2	<code>recipientAddress</code>	-	string	Spinnereistraße 7, 95445 Bayreuth
3	<code>receivedAt</code>	-	string	2022-11-11T16:14:43.692Z

## format

The `format` object specifies the definition or description of the form or characteristics the material should be provided (e.g. to fulfill the order) and represents a(n arbitrary) definition of sample characteristics. The values for the attributes serve as a template from which the resulting material sample(s) in the order typically deviate.

#	Parameter	Parent	Type	Example
1	<code>format</code>	-	<b>object</b>	-
2	<code>id</code>	<code>format</code>	string	-
3	<code>name</code>	<code>format</code>	string	Cypress Format
4	<code>targetContainerType</code>	<code>format</code>	string	GR384
5	<code>targetBarcode</code>	<code>format</code>	string	-
6	<code>targetAmount</code>	<code>format</code>	number	-
7	<code>materialState</code>	<code>format</code>	object	Liquid
8	<code>wellVectors</code>	<code>format</code>	<b>object (array)</b>	?
9	<code>id</code>	<code>wellVectors</code>	string	-
10	<code>instanceNumber</code>	<code>wellVectors</code>	number	0
11	<code>type</code>	<code>wellVectors</code>	string	Sample
12	<code>index</code>	<code>wellVectors</code>	number	?
13	<code>start</code>	<b><code>wellVectors</code></b>	<b>object</b>	-
14	<code>row</code>	<code>start</code>	number	0
15	<code>column</code>	<code>start</code>	number	0
16	<code>wells</code>	<code>wellVectors</code>	number	1
17	<code>direction</code>	<code>wellVectors</code>	string	East
18	<code>spacing</code>	<code>wellVectors</code>	number	0
19	<code>dilutionSeries</code>	<code>wellVectors</code>	array	[ ]
20	<code>rows</code>	<code>wellVectors</code>	number	16
21	<code>columns</code>	<code>wellVectors</code>	number	24

## lineItems

The `lineItems` object breaks a single order down into its component parts and recipients (for notification and/or shipping purposes). `lineItems` determine the parts of an order (expressed as the resulting material sample) that must be processed together and therein attributes.

#	Parameter	Parent	Type	Example
---	-----------	--------	------	---------

1	lineItems	-	<b>object</b>	-
2	id	lineItems	string	-
3	status	lineItems	string	-
4	amount	lineItems	number	0
5	topConcentration	lineItems	number	-
6	recipients	lineItems	array	-
7	recipientAddress	lineItems	string	-

## intermediateFormat

The `intermediateFormat` object is identical to the standard `format` in that it represents the set of characteristics that define the material sample. The difference is that this object enables the material sample final result to be processed in fragments or separate parts.

#	Parameter	Parent	Type	Example
1	<code>intermediateFormat</code>	lineItems	<b>object-format (array)</b>	-
2	id	intermediateFormat	string	-
3	name	intermediateFormat	string	-
4	targetContainertype	intermediateFormat	string	-
5	targetBarcode	intermediateFormat	string	-
6	targetAmount	intermediateFormat	number	-
7	materialState	intermediateFormat	<b>object</b>	-

## comment

The `comment` parameter differs from the `comments` object in that it is a string and is a child parameter of `lineItem`.

#	Parameter	Parent	Type	Example
1	comment	lineItems	string	-
2	layoutIndex	lineItems	number	0
3	deliveryIndex	lineItems	number	0
4	externalReferences	lineItems	string (array)	-
5	sourceType	lineItems	<b>object</b>	-
6	deliverableSampleID	lineItems	string	-

## miracleRun

The `miracleRun` object specifies sample management processing using an intermediate format template.

#	Parameter	Parent	Type	Example
1	miracleRun	-	<b>object</b>	-
2	id	miracleRun	string	-
3	orderID	miracleRun	string	-
4	format	miracleRun	<b>object</b>	-
5	wellVectors	format	<b>object</b>	-
6	amount	miracleRun	string	0
7	topConcentration	miracleRun	string	0
8	equipments	miracleRun	array[string]	-
9	instructions	miracleRun	string	-

## Nested Attributes

The following objects and parameters specify sample management logic specific to Lab Operations.

### wellVectors

The wellVectors object is a `format` attribute that specifies the starting point, type, direction, spacing and an index of the format layout.

✓ Can be a child of `format`, `IntermediateFormat`

#	Parameter	Parent	Type	Example
1	wellVectors	intermediateFormat	object (array)	-
2	id	wellVectors	string	-
3	instanceNumber	wellVectors	number	-
4	type	wellVectors	object	-
5	index	wellVectors	number	-
6	start	wellVectors	object	-
7	row	start	number	-
8	column	start	number	-
9	wells	wellVectors	number	-
10	direction	wellVectors	object	-
11	spacing	wellVectors	number	-
12	dilutionSeries	intermediateFormat	array	[ null ]
13	rows	intermediateFormat	number	16
14	columns	intermediateFormat	number	24

### sourceInventoryItems

The sourceInventoryItems object identifies the material sample in the lab inventory, including the attributes `containerBarcode`, `displayName`, `capacity`, `hazard flags`, `parentID` and `children` among others.

✓ Child of `lineItems`

#	Parameter	Parent	Type	Example
1	sourceInventoryItem	lineItems	object	-
2	sampleID	sourceInventoryItem	string	-
3	containerBarcode	sourceInventoryItem	string	-
4	location	sourceInventoryItem	object	-
5	id	location	string	-
6	displayName	location	string	-
7	type	location	object	-
8	parentID	location	string	-
9	barcode	location	string	-
10	capacity	location	number	0
11	layout	location	array[string]	-
12	orientation	location	object	-
13	counting direction	location	object	-
14	loadingSequence	location	object	-
15	children	location	array	-
16	locationTree	location	array	-
17	hidden	location	boolean	true

18	allowHazardous	location	boolean	true
19	allowControlled	location	boolean	true
20	position	sourceInventoryItem	<b>object</b>	-
21	row	sourceInventoryItem	number	0
22	column	sourceInventoryItem	number	0

## orderFragments

The `orderFragmentation` object is an order segment as a job item. This job item is processed using the rummage table.

<sup>3</sup> The parameter `isLineItemSequenceRelevant` is a boolean data type that is typically nested within an order as an attribute of an `orderFragment`. It can also exist outside of the `orderFragment` parameter.

▼ [Can be a child of lineItems](#)

#	Parameter	Parent	Type	Example
1	orderFragments	-	<b>object (array)</b>	-
2	id	orderFragments	string	-
3	order	<b>orderFragments</b>	<b>object</b>	-
4	id	order	string	-
5	name	order	string	-
6	type	order	object	-
7	description	order	string	-
8	recipients	order	array[string]	-
9	priority	order	object	-
10	receivedAt	order	date	2022-11-21T12:09:44.800Z
11	recipientAddress	order	string	-
12	numberOfCopies	order	object	-
13	<sup>3</sup> isLineItemSequenceRelevant	order	boolean	true

## selectedPath

The `selectedPath` object specifies the processing path (equipment, process, format and job items) selected for the sample management processing sequence.

▼ [Child of orderFragments](#)

#	Parameter	Parent	Type	Example
1	selectedPath	orderFragments	<b>object</b>	-
2	consumed	selectedPath	number	0
3	effort	selectedPath	number	0
4	runs	selectedPath	array[string]	-
5	strategy	selectedPath	object	-